

Product datasheet for **GA200960**

Mouse *Cx3cr1* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cx3cr1 CRISPRa kit - CRISPR gene activation of mouse chemokine (C-X3-C motif) receptor 1
Symbol:	Cx3cr1
UniProt ID:	<u>Q9Z0D9</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)