

## Product datasheet for GA200976

### Mouse Cyp2a12 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| <b>Description:</b>  | Cyp2a12 CRISPRa kit – CRISPR gene activation of mouse cytochrome P450, family 2, subfamily a, polypeptide 12                                                                                                                                                                                                     |
| <b>Symbol:</b>       | Cyp2a12                                                                                                                                                                                                                                                                                                          |
| <b>Components:</b>   | GA200976G1, Cyp2a12 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>GA200976G2, Cyp2a12 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>GA200976G3, Cyp2a12 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>1 CRISPRa-Enhancer vector, SKU GEI00056<br>1 CRISPRa scramble vector, SKU GEI00077 |
| <b>UniProt ID:</b>   | <a href="#">P56593</a>                                                                                                                                                                                                                                                                                           |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:

