

## Product datasheet for GA200985

### Mouse Cyp2c38 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| <b>Description:</b>  | Cyp2c38 CRISPRa kit - CRISPR gene activation of mouse cytochrome P450, family 2, subfamily c, polypeptide 38                                                                                                                                                                                                     |
| <b>Symbol:</b>       | Cyp2c38                                                                                                                                                                                                                                                                                                          |
| <b>Components:</b>   | GA200985G1, Cyp2c38 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA200985G2, Cyp2c38 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA200985G3, Cyp2c38 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <a href="#">P56655</a>                                                                                                                                                                                                                                                                                           |
| <b>Synonyms:</b>     | MGC183632                                                                                                                                                                                                                                                                                                        |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:

