

Product datasheet for **GA200986**

Mouse Cyp2c39 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cyp2c39 CRISPRa kit - CRISPR gene activation of mouse cytochrome P450, family 2, subfamily c, polypeptide 39
Symbol:	Cyp2c39
Components:	GA200986G1, Cyp2c39 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA200986G2, Cyp2c39 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA200986G3, Cyp2c39 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P56656</u>
Synonyms:	CYP11C39
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

