

Product datasheet for **GA200987**

Mouse Cyp2c40 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cyp2c40 CRISPRa kit - CRISPR gene activation of mouse cytochrome P450, family 2, subfamily c, polypeptide 40
Symbol:	Cyp2c40
Components:	GA200987G1, Cyp2c40 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA200987G2, Cyp2c40 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA200987G3, Cyp2c40 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>P56657</u>
Synonyms:	AI265721; AI662255
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

