

Product datasheet for GA201000

Mouse Cyp4a14 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cyp4a14 CRISPRa kit – CRISPR gene activation of mouse cytochrome P450, family 4, subfamily a, polypeptide 14
Symbol:	Cyp4a14
Components:	GA201000G1, Cyp4a14 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA201000G2, Cyp4a14 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA201000G3, Cyp4a14 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q35728</u>
Synonyms:	AI314743
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

