

Product datasheet for **GA201016**

Mouse Dapk3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dapk3 CRISPRa kit - CRISPR gene activation of mouse death-associated protein kinase 3
Symbol:	Dapk3
UniProt ID:	<u>O54784</u>
Synonyms:	dlk; ZIPK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)