

Product datasheet for **GA201085**

Mouse Dlx3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dlx3 CRISPRa kit - CRISPR gene activation of mouse distal-less homeobox 3
Symbol:	Dlx3
UniProt ID:	Q64205
Synonyms:	AV237891; Dlx-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)