

Product datasheet for **GA201086**

Mouse Dlx4 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Dlx4 CRISPRa kit - CRISPR gene activation of mouse distal-less homeobox 4 |
| Symbol: | Dlx4 |
| UniProt ID: | P70436 |
| Synonyms: | Dlx-4; Dlx7 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)