

Product datasheet for **GA201103**

Mouse Dync1i1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dync1i1 CRISPRa kit - CRISPR gene activation of mouse dynein cytoplasmic 1 intermediate chain 1
Symbol:	Dync1i1
Components:	GA201103G1, Dync1i1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA201103G2, Dync1i1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA201103G3, Dync1i1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q88485</u>
Synonyms:	DH IC-1; DHIC-1; DIC; Dncic1; IC74
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)