

Product datasheet for **GA201131**

Mouse Arid3a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Arid3a CRISPRa kit - CRISPR gene activation of mouse AT rich interactive domain 3A (BRIGHT-like)
Symbol:	Arid3a
UniProt ID:	Q62431
Synonyms:	Bright; Dri1; Dril1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)