

Product datasheet for GA201150

Mouse Usp17la activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Usp17la CRISPRa kit - CRISPR gene activation of mouse ubiquitin specific peptidase 17-like A
Symbol:	Usp17la
Components:	GA201150G1, Usp17la gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA201150G2, Usp17la gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA201150G3, Usp17la gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q61068
Synonyms:	Dub1; USP17-A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

