

Product datasheet for GA201151

Mouse Usp17lc activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Usp17lc CRISPRa kit – CRISPR gene activation of mouse ubiquitin specific peptidase 17-like C
Symbol:	Usp17lc
Components:	GA201151G1, Usp17lc gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA201151G2, Usp17lc gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA201151G3, Usp17lc gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>G5E8I7</u>
Synonyms:	Dub-2; Dub2; Dub2b; Usp17I5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

