

Product datasheet for GA201214

Mouse Eif4a1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Eif4a1 CRISPRa kit - CRISPR gene activation of mouse eukaryotic translation initiation factor 4A1
Symbol:	Eif4a1
UniProt ID:	P60843
Synonyms:	BM-010; Ddx2a; Eif4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)