

Product datasheet for **GA201253**

Mouse Epb41l4a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Epb41l4a CRISPRa kit - CRISPR gene activation of mouse erythrocyte membrane protein band 4.1 like 4a
Symbol:	Epb41l4a
Components:	GA201253G1, Epb41l4a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA201253G2, Epb41l4a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA201253G3, Epb41l4a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P52963</u>
Synonyms:	Epb4.1l4; Epb4.1l4a; NBL4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)