

Product datasheet for **GA201357**

Mouse FasI activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FasI CRISPRa kit - CRISPR gene activation of mouse Fas ligand (TNF superfamily, member 6)
Symbol:	FasI
UniProt ID:	<u>P41047</u>
Synonyms:	APT1LG1; CD95-L; CD95L; CD178; Fas-L; Faslg; gld; Tnfsf6; Tnlg1a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)