

Product datasheet for **GA201443**

Mouse Fli1 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Fli1 CRISPRa kit - CRISPR gene activation of mouse Friend leukemia integration 1 |
| Symbol: | Fli1 |
| UniProt ID: | P26323 |
| Synonyms: | EWSR2; Fli-1; SIC-1; Sic1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)