

Product datasheet for **GA201540**

Mouse Galr3 activation kit by CRISPRa

Product data:

| | |
|---------------|-------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | Galr3 CRISPRa kit - CRISPR gene activation of mouse galanin receptor 3 |
| Symbol: | Galr3 |
| UniProt ID: | <u>O88853</u> |
| Synonyms: | Galnr3 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)