

## Product datasheet for **GA201577**

### Mouse Mrps33 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Mrps33 CRISPRa kit - CRISPR gene activation of mouse mitochondrial ribosomal protein S33
Symbol:	Mrps33
Components:	GA201577G1, Mrps33 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA201577G2, Mrps33 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA201577G3, Mrps33 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D2R8</u>
Synonyms:	AI841153; CGI-139; Gdap3; MRP-S33; PTD003
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)