

Product datasheet for **GA201620**

Mouse **Gjb3** activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Gjb3 CRISPRa kit - CRISPR gene activation of mouse gap junction protein, beta 3 |
| Symbol: | Gjb3 |
| UniProt ID: | P28231 |
| Synonyms: | Cnx31; Cx31; D4Wsu144e; Gjb-3 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)