

Product datasheet for **GA201644**

Mouse **Gnl1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Gnl1 CRISPRa kit - CRISPR gene activation of mouse guanine nucleotide binding protein-like 1
Symbol:	Gnl1
Components:	GA201644G1, Gnl1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA201644G2, Gnl1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA201644G3, Gnl1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	P36916
Synonyms:	Gna-rs1; Gnal1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)