

Product datasheet for **GA201668**

Mouse Gng3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Gng3 CRISPRa kit - CRISPR gene activation of mouse guanine nucleotide binding protein (G protein), gamma 3
Symbol:	Gng3
UniProt ID:	P63216
Synonyms:	83
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)