

Product datasheet for **GA201686**

Mouse **Lilrb4a** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Lilrb4a CRISPRa kit - CRISPR gene activation of mouse leukocyte immunoglobulin-like receptor, subfamily B, member 4A
Symbol:	Lilrb4a
UniProt ID:	Q64281
Synonyms:	CD85K; gp49; Gp49b; HM18; ILT3; Lilrb4; LIR-5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)