

Product datasheet for **GA201704**

Mouse Gpr33 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Gpr33 CRISPRa kit - CRISPR gene activation of mouse G protein-coupled receptor 33
Symbol:	Gpr33
Components:	GA201704G1, Gpr33 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA201704G2, Gpr33 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA201704G3, Gpr33 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O88416</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)