

## Product datasheet for GA201725

### Mouse Grcc10 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | Grcc10 CRISPRa kit – CRISPR gene activation of mouse gene rich cluster, C10 gene                                                                                                                                                                                                                              |
| <b>Symbol:</b>       | Grcc10                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA201725G1, Grcc10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>GA201725G2, Grcc10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>GA201725G3, Grcc10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>1 CRISPRa-Enhancer vector, SKU GEI00056<br>1 CRISPRa scramble vector, SKU GEI00077 |
| <b>UniProt ID:</b>   | <u><a href="#">O35127</a></u>                                                                                                                                                                                                                                                                                 |
| <b>Synonyms:</b>     | 2310033H05Rik; C10                                                                                                                                                                                                                                                                                            |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

