

Product datasheet for GA201750

Mouse Pdia3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pdia3 CRISPRa kit - CRISPR gene activation of mouse protein disulfide isomerase associated 3
Symbol:	Pdia3
UniProt ID:	P27773
Synonyms:	58kDa; Erp; ERp57; ERp60; ERp61; Grp58; PDI; PDI-Q2; PI-PLC; Plca; PLC[a]
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)