

Product datasheet for **GA201832**

Mouse *Ifi44l* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ifi44l CRISPRa kit - CRISPR gene activation of mouse interferon-induced protein 44 like
Symbol:	Ifi44l
UniProt ID:	Q9BDB7
Synonyms:	H-28; H28; H28-1; NS1178
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)