

Product datasheet for **GA201836**

Mouse H3f3b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	H3f3b CRISPRa kit - CRISPR gene activation of mouse H3 histone, family 3B
Symbol:	H3f3b
UniProt ID:	P84244
Synonyms:	9430068D06Rik; H3.3B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)