

Product datasheet for **GA201877**

Mouse Hes3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hes3 CRISPRa kit - CRISPR gene activation of mouse hes family bHLH transcription factor 3
Symbol:	Hes3
UniProt ID:	Q61657
Synonyms:	bHLHb43; Hes-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)