

Product datasheet for **GA201903**

Mouse Hipk3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hipk3 CRISPRa kit - CRISPR gene activation of mouse homeodomain interacting protein kinase 3
Symbol:	Hipk3
UniProt ID:	Q9ERH7
Synonyms:	DYRK6; FIST3; mir-1902; PKY
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)