

Product datasheet for GA201970

Mouse Hoxd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hoxd3 CRISPRa kit - CRISPR gene activation of mouse homeobox D3
Symbol:	Hoxd3
Components:	GA201970G1, Hoxd3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA201970G2, Hoxd3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA201970G3, Hoxd3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P09027</u>
Synonyms:	Hox-4.1; Hox-5.5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

