

Product datasheet for GA201975

Mouse Hp1bp3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hp1bp3 CRISPRa kit - CRISPR gene activation of mouse heterochromatin protein 1, binding protein 3
Symbol:	Hp1bp3
UniProt ID:	Q3TEA8
Synonyms:	Hp1bp74
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)