

Product datasheet for **GA202068**

Mouse *Ifi47* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ifi47 CRISPRa kit - CRISPR gene activation of mouse interferon gamma inducible protein 47
Symbol:	Ifi47
Synonyms:	47kDa; Ifggc1; Igrd; ligp4; lipg4; IRG-47; Irgd
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)