

Product datasheet for **GA202069**

Mouse Ifit2 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Ifit2 CRISPRa kit - CRISPR gene activation of mouse interferon-induced protein with tetratricopeptide repeats 2 |
| Symbol: | Ifit2 |
| UniProt ID: | Q64112 |
| Synonyms: | AV302338; Ifi54; P54 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)