

Product datasheet for **GA202078**

Mouse *Ifnar1* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ifnar1 CRISPRa kit - CRISPR gene activation of mouse interferon (alpha and beta) receptor 1
Symbol:	Ifnar1
UniProt ID:	P33896
Synonyms:	CD118; Ifar; Ifnar; Ifrc; Infar
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)