

Product datasheet for **GA202109**

Mouse Cd74 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cd74 CRISPRa kit - CRISPR gene activation of mouse CD74 antigen (invariant polypeptide of major histocompatibility complex, class II antigen-associated)
Symbol:	Cd74
UniProt ID:	P04441
Synonyms:	CLIP; DHLAG; HLADG; Ia-GAMMA; li
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)