

Product datasheet for **GA202129**

Mouse *Il17ra* activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Il17ra CRISPRa kit - CRISPR gene activation of mouse interleukin 17 receptor A |
| Symbol: | Il17ra |
| UniProt ID: | Q60943 |
| Synonyms: | AW538159; Cdw217; Il17r; VDw217 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)