

Product datasheet for **GA202211**

Mouse Cd47 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cd47 CRISPRa kit - CRISPR gene activation of mouse CD47 antigen (Rh-related antigen, integrin-associated signal transducer)
Symbol:	Cd47
UniProt ID:	Q61735
Synonyms:	9130415E20Rik; AA407862; AI848868; AW108519; B430305P08Rik; IAP; Itgp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)