

Product datasheet for **GA202215**

Mouse *Itih4* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Itih4 CRISPRa kit - CRISPR gene activation of mouse inter alpha-trypsin inhibitor, heavy chain 4
Symbol:	Itih4
Components:	GA202215G1, Itih4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA202215G2, Itih4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA202215G3, Itih4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A6X935
Synonyms:	ITI-HC4; Itih-4; PK-120
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)