

Product datasheet for **GA202243**

Mouse Kap activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kap CRISPRa kit - CRISPR gene activation of mouse kidney androgen regulated protein
Symbol:	Kap
Components:	GA202243G1, Kap gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA202243G2, Kap gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA202243G3, Kap gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P61110</u>
Synonyms:	AW146435; FKSG22
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)