

Product datasheet for **GA202244**

Mouse *Kcna1* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcna1</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, shaker-related subfamily, member 1
Symbol:	<i>Kcna1</i>
UniProt ID:	P16388
Synonyms:	AI840627; Kca1-1; Kv1.1; MBK1; mceph; Mk-1; Shak
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)