

Product datasheet for **GA202245**

Mouse *Kcna2* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcna2</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, shaker-related subfamily, member 2
Symbol:	<i>Kcna2</i>
UniProt ID:	P63141
Synonyms:	Akr6a4; ENSMUSG00000074335; Gm10672; Kca1-2; Kv1.2; Mk-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)