

Product datasheet for **GA202246**

Mouse *Kcna3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcna3</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, shaker-related subfamily, member 3
Symbol:	<i>Kcna3</i>
UniProt ID:	P16390
Synonyms:	Kca1-3; Kv1.3; Mk-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)