

## Product datasheet for **GA202249**

### Mouse **Kcna6** activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Kcna6 CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, shaker-related, subfamily, member 6
Symbol:	Kcna6
UniProt ID:	<a href="#">Q61923</a>
Synonyms:	Kv1.6; MK1.6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)