

Product datasheet for **GA202251**

Mouse *Kcnab1* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnab1</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, shaker-related subfamily, beta member 1
Symbol:	<i>Kcnab1</i>
UniProt ID:	P63143
Synonyms:	Akr8a8; Kvbeta1.1; mKv(beta)1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)