

Product datasheet for **GA202253**

Mouse **Kcnab3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Kcnab3 CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, shaker-related subfamily, beta member 3
Symbol:	Kcnab3
Synonyms:	C330022D06Rik; Kcnab4; mKv(beta)4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)