

Product datasheet for **GA202254**

Mouse **Kcnb1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Kcnb1 CRISPRa kit - CRISPR gene activation of mouse potassium voltage gated channel, Shab-related subfamily, member 1
Symbol:	Kcnb1
UniProt ID:	Q03717
Synonyms:	Kcr1-1; Kv2.1; Shab
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)