

## Product datasheet for **GA202255**

### Mouse *Kcnc1* activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnc1</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage gated channel, Shaw-related subfamily, member 1
Symbol:	<i>Kcnc1</i>
UniProt ID:	<a href="#">P15388</a>
Synonyms:	C230009H10Rik; <i>Kcr2-1</i> ; <i>KShIIIB</i> ; <i>Kv3.1</i> ; <i>KV4</i> ; <i>NGK2</i> ; Shaw
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)