

Product datasheet for **GA202256**

Mouse *Kcnc3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnc3</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage gated channel, Shaw-related subfamily, member 3
Symbol:	<i>Kcnc3</i>
UniProt ID:	Q63959
Synonyms:	Kcr2-3; KShIIID; Kv3.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)