

## Product datasheet for **GA202260**

### Mouse *Kcnh1* activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnh1</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, subfamily H (eag-related), member 1
Symbol:	<i>Kcnh1</i>
UniProt ID:	<a href="https://www.uniprot.org/entry/Q60603">Q60603</a>
Synonyms:	EAG1; Kv10.1; M-eag
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)